



Product description

| | |
|---------------------|---|
| Order number | 75.570 |
| Product description | Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, Ø opening: 70 mm, white, cap: transparent, material: PS, snap-on lid, 50 piece(s)/tube bag |

Product characteristics

| | |
|----------------|-------------|
| Graduation | no |
| Closure type | snap-on lid |
| Volume of work | 100 ml |

Size

| | |
|----------------------|--------|
| Sample volume | 100 ml |
| Max. volume | 100 ml |
| Diameter | 70 mm |
| Diameter of opening | 70 mm |
| Height of product | 49 mm |
| Length of product | 49 mm |
| Length excluding cap | 49 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, PS, white

Page 2

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | white |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | tube bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 45000 |
| Case volume | 0.0317 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 3.71 kg |
| EAN of inner box | 4038917082563 |
| EAN case | 4038917001649 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-22

This document was prepared by EDP support and is valid without signature.



SARSTEDT