



## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, translucent,  
with paper label**

**Page 1**



### Product description

Order number	75.9922.730
Product description	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, translucent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: yellow, sterile, 250 piece(s)/bag

### Product characteristics

Label/ Print	with paper label
Colour of print/label	yellow
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	50 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, translucent,  
with paper label**

**Page 2**

### Size

Max. volume	70 ml
Diameter	44 mm
Height of product	55 mm
Length of product	55 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	15000
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of case	7 kg
EAN of inner box	4038917152617
EAN case	4038917038836

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**