



Product description

Order number	76.1477.100
Product description	Snap-on lid, with spout, suitable for 125-ml container 75.1355, material: PE, red, 100 piece(s)/bag

Product characteristics

Closure type	snap-on lid
Design	with spout
Suitable for	125-ml container 75.1355

Material & colours

Product material	Polyethylene (PE)
Colour of product	red

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	59 mm
Width of case	40 mm
Height of case	28 mm
Case volume	0.0646 cbm
Weight of product	0.0037 kg
Weight of case	4.57 kg
EAN of inner box	4038917129411
EAN case	4042836001604

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

