



Product description

| | |
|---------------------|--|
| Order number | 76.564 |
| Product description | Screw cap, suitable for 100-ml container 75.563, material: HD-PE, natural, 50 piece(s)/bag |

Product characteristics

| | |
|----------------|-------------------------------|
| Closure type | screw cap |
| Example of use | cap for urine collection tube |
| Suitable for | 100-ml container 75.563 |

Size

| | |
|-------------------|-------|
| Diameter | 62 mm |
| Height of product | 14 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of product | natural |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 21000 |
| Depth of case | 388 mm |
| Width of case | 378 mm |
| Height of case | 250 mm |
| Case volume | 0.0367 cbm |
| Weight of product | 0.0068 kg |
| Weight of case | 3.91 kg |
| EAN of inner box | 4038917082518 |
| EAN case | 4038917001595 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

