



Product description

| | |
|---------------------|---|
| Order number | 51.9923.820 |
| Product description | Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, skirted conical base, transparent, cap assembled, with paper label, preparation: Boric acid, non-sterile, 500 piece(s)/case |

Product characteristics

| | |
|--------------------------|----------------------|
| Type of preparation | Boric acid |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 25 ml |
| Centrifugation max (RCF) | 4500 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile

Page 2

Size

| | |
|-------------------|-------|
| Sample volume | 25 ml |
| Diameter | 25 mm |
| Length of product | 90 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Colour of cap | red |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | non-sterile |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 15000 |
| Case volume | 0.059 cbm |
| Weight of product | 0.009 kg |
| Weight of case | 5.18 kg |
| EAN of inner box | 4038917154550 |
| EAN case | 4038917030090 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT