



### **Product description**

Order number 57.527  
Product description Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PP, conical base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Base shape conical base  
Volume of work 12 ml  
Centrifugation max (RCF) 3000 x g

### **Size**

Diameter 17 mm  
Length of product 110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0052 kg
Weight of case	5.93 kg
EAN of inner box	4038917082228
EAN case	4038917001311

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

