



### **Product description**

Order number	60.9921.829
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PS, skirted conical base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	skirted conical base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16 mm
Length of product	97 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### **Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0077 kg
Weight of case	8.5 kg
EAN of inner box	4038917151122
EAN case	4038917027700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

