




Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**



Product description

| | |
|---------------------|---|
| Order number | 62.554.502 |
| Product description | Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Label/ Print | with print |
| Colour of print/label | white/blue |
| Graduation | yes |
| Cap | standard screw cap |
| Closure type | screw cap |
| Base shape | conical base |
| Volume of work | 15 ml |
| Centrifugation max (RCF) | 20000 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Diameter | 17 mm |
| Length of product | 120 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | red |

Purity & certification

| | |
|---------------------------|--|
| Satisfies the requirement | IATA |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 20000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |
| Height of case | 165 mm |
| Case volume | 0.0393 cbm |
| Weight of product | 0.0063 kg |
| Weight of case | 3.74 kg |
| EAN of inner box | 4038917113182 |
| EAN case | 4038917013932 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

