



Product description

Order number	10.251
Product description	Collection tip for Urine Monovette® , for Urine transfer into Monovette® Urine, 100 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Light protection	no
Application	Urine transfer into Monovette® Urine
Suitable for	Urine Monovette®
Connection	Luer (f)

Size

Diameter	8 mm
Length of product	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product yellow

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 5000
Piece(s) / pallet 240000
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0005 kg
Weight of case 3.31 kg
EAN of inner box 4038917080941
EAN case 4038917000116

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

