



Product Specification

Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 1



Product description

Order number	10.253.388
Product description	Urine Monovette® , Boric acid, nominal volume: 10 ml, Boric acid, Urine tube, (LxØ): 102 x 15 mm, transparent, material: PP, Luer screw cap, cap green, incl. collection tip, with plastic label pre-barcode, Pre-barcode with 8-digit unique number range and 3-digit prefix, sterile, 1 piece(s)/blister

Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with plastic label pre-barcode
Colour of print/label	white/green
Graduation	no
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

Page 2

Volume of work 10 ml
Connection Luer (m)

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 117 mm
Length excluding cap 102 mm
Length including cap and plunger 127 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 20000
Type of smallest blister
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 500
Piece(s) / pallet 13500
Depth of box 351 mm
Width of box 151 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 1 piece(s)/blister**

Page 3

Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0085 kg
Weight of case	5.89 kg
EAN of inner box	4038917950008
EAN case	4038917950015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT