



### Product description

Order number	10.258.001
Product description	Urine Monovette® , Z, nominal volume: 8.5 ml, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, 64 piece(s)/bag

### Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow/white
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	8,5 ml
Recommended centrifugation	400 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Z, 8.5 ml, cap yellow, (LxØ): 92  
x 15 mm, 64 piece(s)/bag**

**Page 2**

Centrifugation time 5 Min  
Connection Luer (m)

**Size**

Sample volume 8.5 ml  
Diameter 15 mm  
Length including cap 106 mm  
Length excluding cap 92 mm  
Length including cap and plunger 127 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 512  
Type of smallest subpackaging bag  
Piece(s) / inner box 64  
Piece (s) / outer case 512  
Piece(s) / pallet 24576  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 215 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Urine Monovette® , Z, 8.5 ml, cap yellow, (LxØ): 92  
x 15 mm, 64 piece(s)/bag**

**Page 3**

Case volume	0.0315 cbm
Weight of product	0.0078 kg
Weight of case	4.49 kg
EAN of inner box	4038917280020
EAN case	4038917279888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**