



Product description

Order number	10.258.061
Product description	Urine Monovette® , Z, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with paper label, suitable for Tempus600, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow/white
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus600
Base shape	round base
Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

Size

Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 64 piece(s)/bag

Page 2

Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 512
Type of smallest bag
subpackaging
Piece(s) / inner box 64
Piece (s) / outer case 512
Piece(s) / pallet 24576
Case volume 0.0315 cbm
Weight of product 0.0078 kg
Weight of case 4.49 kg
EAN of inner box 4038917992435
EAN case 4038917992442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT