



## Product description

Order number	10.258.081
Product description	Urine Monovette® , Z, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, Luer screw cap, cap yellow, incl. collection tip, with plastic label, suitable for Tempus600, 1 piece(s)/blister

## Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label
Colour of print/label	yellow/white
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus600
Base shape	round base
Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

## Size

Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

## Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Z, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 2**

Piston plunger material                      Polystyrene (PS)

**Purity & certification**

Satisfies the requirement                      IATA, ADR  
Product category                                      In vitrodiagnostic device, CE  
CE certified    CE  
Sterilisation    Electron irradiation  
Safety data sheet                                      no  
Batched    yes

**Packaging**

Minimum order qty.                                      500  
Type of smallest  
subpackaging    blister  
Piece(s) / inner box                                      100  
Piece (s) / outer case                                      500  
Piece(s) / pallet    13500  
Depth of box    351 mm  
Width of box    151 mm  
Height of box    171 mm  
Depth of case    778 mm  
Width of case    368 mm  
Height of case    195 mm  
Case volume    0.0558 cbm  
Weight of product    0.0078 kg  
Weight of case    5.51 kg  
EAN of inner box    4038917992558  
EAN case    4038917992565

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-08-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**