




**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
orange, with handle**

**Page 1**



**Product description**

Order number	93.851.175
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, orange, with handle, material: PP

**Product characteristics**

Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

**Size**

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	246 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
orange, with handle**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product orange

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 12  
Piece(s) / pallet 840  
Depth of case 453 mm  
Width of case 258 mm  
Height of case 165 mm  
Case volume 0.0193 cbm  
Weight of product 0.1983 kg  
Weight of case 2.79 kg  
EAN of inner box 4038917121873  
EAN case 4038917075015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**