



**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
creamy white**

**Page 1**



### **Product description**

Order number	93.852.170
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, creamy white, material: recycled PP

### **Product characteristics**

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

### **Size**

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Block Rack D17, Ø opening: 17 mm, 5 x 10,  
creamy white**

**Page 2**

**Material & colours**

Product material recycled Polypropylene (recycled PP)  
Colour of product creamy white

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 760  
Depth of case 568 mm  
Width of case 433 mm  
Height of case 120 mm  
Case volume 0.0295 cbm  
Weight of product 0.1958 kg  
Weight of case 4.62 kg  
EAN of inner box 4038917122054  
EAN case 4038917044905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**