



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 93.852.300   |
| Product description | Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 5 x 20, grey, material: recycled PP |

### **Product characteristics**

|                        |        |
|------------------------|--------|
| Number of wells        | 100    |
| Storage capacity       | 100    |
| Format                 | 5 x 20 |
| Autoclavable (max. °C) | 121 °C |

### **Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 17 mm  |
| Width of product    | 109 mm |
| Height of product   | 45 mm  |
| Length of product   | 418 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material recycled Polypropylene (recycled PP)  
Colour of product grey

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 12  
Piece(s) / pallet 528  
Depth of case 453 mm  
Width of case 433 mm  
Height of case 165 mm  
Case volume 0.0324 cbm  
Weight of product 0.3917 kg  
Weight of case 5.41 kg  
EAN of inner box 4038917122221  
EAN case 4038917018296

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

