



### **Product description**

Order number	73.1053
Product description	Insert, suitable for tubes and S-Monovette® Ø 16 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Suitable for	tubes and S-Monovette® Ø 16 mm
Base shape	skirted conical base

### **Size**

Diameter of opening	11 mm
Length of product	23.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Insert, suitable for tubes and S-Monovette® Ø 16 mm, transparent**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 2000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 156000  
Depth of case 428 mm  
Width of case 323 mm  
Height of case 130 mm  
Case volume 0.018 cbm  
Weight of product 0.0009 kg  
Weight of case 2.26 kg  
EAN of inner box 4038917085168  
EAN case 4038917003841

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**