



### **Product description**

Order number	73.667
Product description	Sample tube, suitable for Abbott analyser AXSYM, material: PS, 2.5 ml, (LxØ): 53 x 12 mm, transparent, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Suitable for	Abbott analyser AXSYM
Base shape	conical base

### **Size**

Max. volume	2.5 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	53 mm
Length of product	53 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 2500  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 2500  
Piece(s) / pallet 90000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 265 mm  
Case volume 0.042 cbm  
Weight of product 0.0025 kg  
Weight of case 6.78 kg  
EAN of inner box 4038917083232  
EAN case 4038917002189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

