




Product Specification
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,
with paper label, transparent**

Page 1



Product description

Order number	75.1435.004
Product description	Urine container, max. working volume: 90 ml, (ØxH): 60 x 65 mm, Ø opening: 65 mm, transparent, graduated, material: PP, screw cap, with paper label, cap assembled, sterile, 400 piece(s)/case

Product characteristics

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Volume of work	90 ml

Size

Sample volume	90 ml
Max. volume	90 ml
Diameter	60 mm
Diameter of opening	65 mm
Height of product	65 mm
Length excluding cap	60 µm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,
with paper label, transparent**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 400
Type of smallest case
subpackaging
Piece(s) / inner box 400
Piece (s) / outer case 400
Piece(s) / pallet 7200
Case volume 0.1062 cbm
Weight of product 0.0201 kg
Weight of case 8.23 kg
EAN of inner box 4038917385220
EAN case 4038917385213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT