



Product Specification
S-Monovette® , 9 ml, for Formalin system, (LxØ):
92 x 16 mm, with paper label

Page 1



Product description

Order number	09.1704.001
Product description	S-Monovette® , up to 9 ml, for Formalin system, (LxØ): 92 x 16 mm, with paper label, 50 piece(s)/bag

Product characteristics

Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Closure type	screw cap
Application	Tissue fixation
Suitable for	Formalin system
Volume of work	9 ml
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® , 9 ml, for Formalin system, (LxØ):
92 x 16 mm, with paper label

Page 2

Size

Max. volume	9 ml
Diameter	16 mm
Length of product	92 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0078 kg
Weight of case	4.39 kg
EAN of inner box	4038917335317
EAN case	4038917335300

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

