



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
cap enclosed, 1,000 piece(s)/bag**

**Page 1**



### Product description

Order number	62.612.009
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with print, label/print: white, with graduations, cap enclosed, natural, 1,000 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Base shape	conical false bottom, rounded tube bottom
Volume of work	5 ml
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
cap enclosed, 1,000 piece(s)/bag**

**Page 2**

### Size

Diameter	15.3 mm
Height of product	92 mm
Length of product	92 mm
Length excluding cap	92 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0492 cbm
Weight of product	0.0054 kg
Weight of case	5.9 kg
EAN of inner box	4038917114929
EAN case	4038917072298

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**