



**Product Specification**  
**Thrombo Plus® , 2 ml, cap red, (LxØ): 75 x 11.5 mm, with print**



**Product description**

Order number	51.334
Product description	Thrombo Plus® , for direct platelet counting, nominal volume: 2 ml, cap red, preparation: Hg salt, push cap, (LxØ): 75 x 11.5 mm, round base, with print, 50 piece(s)/case

**Product characteristics**

Type of collection	Venous and capillary
Type of preparation	Hg salt
Preparation concentration	2.5 mg/ml solution
Label/ Print	with print
Colour of print/label	red
Graduation	yes
Closure type	push cap
Application	venous and capillary blood samples
Example of use	for direct platelet counting
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Thrombo Plus® , 2 ml, cap red, (LxØ): 75 x 11.5 mm, with print**

**Page 2**

**Size**

Sample volume	2 ml
Diameter	11.5 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	red

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	48000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0055 kg
Weight of case	3.03 kg
EAN of inner box	4038917081207
EAN case	4038917000345

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**