



### **Product description**

Order number	73.834
Product description	Cuvette and steel balls, suitable for Amelung KC coagulometer (macro), material: PS, length: 23 mm, transparent, 1,000 cuvettes/bag and 1,000 steel balls/flask

### **Product characteristics**

Graduation	no
Suitable for	Amelung KC coagulometer (macro)
Base shape	flat base

### **Size**

Diameter of opening	10 mm
Length of product	23 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

**Packaging**

Minimum order qty.                      2000  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      2000  
Piece(s) / pallet                          216000  
Depth of case                              308 mm  
Width of case                                223 mm  
Height of case                               195 mm  
Case volume                                0.0134 cbm  
Weight of product                        0.0008 kg  
Weight of case                              1.8 kg  
EAN of inner box                          4038917084321  
EAN case                                      4038917003186

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

