



### **Product description**

Order number 73.834.001  
Product description Steel ball for Amelung KC coagulometer (macro), grey

### **Product characteristics**

Suitable for Amelung KC coagulometer (macro)

### **Material & colours**

Colour of product grey

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Piece(s) / inner box 1000  
Piece (s) / outer case 50000  
Piece(s) / pallet 5400000  
Depth of case 308 mm  
Width of case 223 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0001 kg
Weight of case	7.15 kg
EAN of inner box	4038917121682
EAN case	4038917017930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

