



### **Product description**

Order number	90.188.055
Product description	Paper roll, for Thermal printer REF. 90.189.720 & 90.1092.720

### **Product characteristics**

Suitable for	Thermal printer REF. 90.189.720 & 90.1092.720
--------------	---

### **Size**

Width of product	57 mm
Length of product	25 m

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	5
Piece (s) / outer case	10
Piece(s) / pallet	3240
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0808 kg
Weight of case	0.91 kg
EAN of inner box	4038917102759
EAN case	4038917008884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

