



Product description

| | |
|---------------------|--|
| Order number | 51.1534.500 |
| Product description | Salivette® Cortisol, with synthetic swab, cap: blue, with paper label, label/print: white/blue, for cortisol determination from saliva, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|------------------------------------|
| Label/ Print | with paper label |
| Colour of print/label | white/blue |
| Closure type | push cap |
| Design | with synthetic swab |
| Application | cortisol determination from saliva |
| Recommended centrifugation | 1000 x g |
| Centrifugation time | 2 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|----------------------|---------|
| Diameter | 16.8 mm |
| Length including cap | 97 mm |
| Length excluding cap | 84 mm |

Material & colours

| | |
|------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | blue |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 28000 |
| Depth of box | 226 mm |
| Width of box | 101 mm |
| Height of box | 119 mm |
| Depth of case | 533 mm |
| Width of case | 263 mm |
| Height of case | 140 mm |
| Case volume | 0.0196 cbm |
| Weight of product | 0.0069 kg |
| Weight of case | 4.09 kg |
| EAN of inner box | 4038917159364 |
| EAN case | 4038917159357 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

