



### **Product description**

Order number	95.978
Product description	Mailing envelope, material: PE, (LxW): 230 x 300 mm, with 50 mm adhesive tape, film thickness: 0.045 mm, 100 piece(s)/case

### **Product characteristics**

Design	with 50 mm adhesive tape
Absorbent liner	no

### **Size**

Width of product	300 mm
Length of product	230 mm

### **Material & colours**

Product material	Polyethylene (PE)
------------------	-------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging case  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 56000  
Depth of case 290 mm  
Width of case 145 mm  
Height of case 400 mm  
Case volume 0.0168 cbm  
Weight of product 0.0057 kg  
Weight of case 6.12 kg  
EAN of inner box 4038917224178  
EAN case 4038917224161

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

