



Product description

| | |
|---------------------|--|
| Order number | 95.1713 |
| Product description | Cross-bottom bag, (LxW): 270 x 270 mm, material: PE, for transport case B 17, 25 piece(s)/case |

Product characteristics

| | |
|--------------|---------------------|
| Suitable for | transport case B 17 |
|--------------|---------------------|

Size

| | |
|-------------------|--------|
| Width of product | 270 mm |
| Length of product | 270 mm |

Material & colours

| | |
|------------------|-------------------|
| Product material | Polyethylene (PE) |
|------------------|-------------------|

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 25 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 25 |
| Piece (s) / outer case | 25 |
| Piece(s) / pallet | 1800 |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0582 kg |
| Weight of case | 1.7 kg |
| EAN of inner box | 4038917087513 |
| EAN case | 4038917065917 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

