



**Product Specification**  
**Mailing container, length: 179 mm, Ø opening: 28 mm, cap enclosed**

**Page 1**



**Product description**

Order number	78.574
Product description	Mailing container, (LxW): 179 x 35 mm, Ø opening: 28 mm, cap enclosed, cap: natural, 50 piece(s)/bag

**Product characteristics**

Cap	round cap
Closure type	screw cap
Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	no
Shape of opening	round
Liquid-tight	yes

**Size**

Diameter of opening	28 mm
Width of product	35 mm
Length of product	179 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Mailing container, length: 179 mm, Ø opening: 28 mm, cap enclosed**

**Page 2**

**Material & colours**

Product material High Density Polyethylene (HD-PE)  
Colour of cap natural

**Purity & certification**

Satisfies the requirement P 650, ADR  
Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 250  
Type of smallest subpackaging bag  
Piece(s) / inner box 50  
Piece (s) / outer case 250  
Piece(s) / pallet 3000  
Depth of case 563 mm  
Width of case 343 mm  
Height of case 500 mm  
Case volume 0.0966 cbm  
Weight of product 0.0201 kg  
Weight of case 5.94 kg  
EAN of inner box 4038917159449  
EAN case 4038917159432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**