



## Product Specification

**Mailing container, transparent, construction: round, with absorbent liner, length: 126 mm, Ø opening: 30 mm, without cap**

Page 1



### Product description

Order number	78.898
Product description	Mailing container, round, transparent, construction: round, with absorbent liner, length: 126 mm, Ø opening: 30 mm, without cap, 50 piece(s)/bag

### Product characteristics

Application	Mailing of samples
Example of use	Secondary packaging according to P 650
Absorbent liner	yes
Shape of opening	round
Liquid-tight	yes

### Size

Diameter of opening	30 mm
Length of product	126 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Mailing container, transparent, construction:  
round, with absorbent liner, length: 126 mm, Ø  
opening: 30 mm, without cap**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	7500
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	0.0145 kg
Weight of case	4.2 kg
EAN of inner box	4038917084543
EAN case	4038917003353

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**