



### **Product description**

|                     |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 71.9921.316                                                                                                                       |
| Product description | Bucket, up to 1 l, (ØxH): 122 x 145 mm, white, cap enclosed, for collection and storage of diagnostic specimens, 24 piece(s)/case |

### **Product characteristics**

|                        |                                                |
|------------------------|------------------------------------------------|
| Application            | collection and storage of diagnostic specimens |
| Autoclavable (max. °C) | 121 °C                                         |

### **Size**

|                   |        |
|-------------------|--------|
| Max. volume       | 1 l    |
| Diameter          | 122 mm |
| Height of product | 145 mm |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | white              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

|                   |     |
|-------------------|-----|
| Latex-free        | yes |
| Pyrogen-free      | no  |
| Safety data sheet | no  |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 24            |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 24            |
| Piece (s) / outer case        | 24            |
| Piece(s) / pallet             | 1152          |
| Depth of case                 | 428 mm        |
| Width of case                 | 358 mm        |
| Height of case                | 234 mm        |
| Case volume                   | 0.0357 cbm    |
| Weight of case                | 2.7 kg        |
| EAN of inner box              | 4038917150620 |
| EAN case                      | 4038917038195 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

