



**Product Specification**  
**Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation**



**Product description**

Order number	53.419
Product description	Seraplas® valve filter, V13, blue, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

**Product characteristics**

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

**Size**

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation**

**Page 2**

**Material & colours**

Product material Polyethylene (PE)  
Colour of product blue

**Purity & certification**

Product category no medical device | IVD  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging case  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 75000  
Depth of box 119 mm  
Width of box 99 mm  
Height of box 98 mm  
Depth of case 513 mm  
Width of case 293 mm  
Height of case 140 mm  
Case volume 0.021 cbm  
Weight of product 0.0016 kg  
Weight of case 2.15 kg  
EAN of inner box 4038917081481  
EAN case 4038917000574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**