



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 55.1577  |
| Product description | Carrier tube, (LxØ): 60 x 11.5 mm, material: PP, transparent, 250 piece(s)/bag |

### **Product characteristics**

|              |             |
|--------------|-------------|
| Closure type | screw cap   |
| Suitable for | Multivette® |
| Base shape   | flat base   |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 11.5 mm |
| Length of product | 60 mm   |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 250  
Piece (s) / outer case 1000  
Piece(s) / pallet 78000  
Depth of case 428 mm  
Width of case 323 mm  
Height of case 130 mm  
Case volume 0.018 cbm  
Weight of product 0.0013 kg  
Weight of case 1.7 kg  
EAN of inner box 4038917087070  
EAN case 4038917005548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

