



### **Product description**

Order number	55.525.003
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, suitable for Microvette® and micro tube with screw cap, 250 piece(s)/bag

### **Product characteristics**

Suitable for	Microvette® and micro tube with screw cap
Base shape	round base
Volume of work	5 ml

### **Size**

Diameter	13 mm
Length of product	75 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      2000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      250  
Piece (s) / outer case                      2000  
Piece(s) / pallet                            72000  
Depth of case                              398 mm  
Width of case                                398 mm  
Height of case                               265 mm  
Case volume                                0.042 cbm  
Weight of product                        0.0026 kg  
Weight of case                              5.7 kg  
EAN of inner box                         4038917109154  
EAN case                                    4038917012362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

