



Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, flat tube bottom, PP, without
cap, 100 piece(s)/bag**

Page 1



Product description

Order number	60.613.010
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, flat tube bottom, transparent, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Graduation	no
Base shape	conical false bottom, flat tube bottom
Volume of work	3,5 ml

Size

Diameter	13 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, flat tube bottom, PP, without
cap, 100 piece(s)/bag**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 42000
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0026 kg
Weight of case 3.37 kg
EAN of inner box 4038917114974
EAN case 4038917079143

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT