



Product description

Order number	65.163
Product description	Screw cap, material: HD-PE, natural, for tubes 75 x 13 mm, 1,000 piece(s)/bag

Product characteristics

Suitable for	tubes 75 x 13 mm
--------------	------------------

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0014 kg
Weight of case	7.45 kg
EAN of inner box	4038917272223
EAN case	4038917272216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

