



Product description

Order number 65.647.320
Product description with filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag

Product characteristics

Suitable for S-Monovette® , tubes Ø 13-16 mm
Filter yes

Material & colours

Product material Polypropylene (PP)
Colour of product light blue

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	378 mm
Width of case	250 mm
Height of case	388 mm
Case volume	0.0367 cbm
Weight of case	5.37 kg
EAN of inner box	4038917988278
EAN case	4038917988285

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

