



Product description

Order number	65.806
Product description	Push cap, natural, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 13 mm

Size

Diameter	13.3 mm
Height of product	15.2 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 10000
Type of smallest
subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 10000
Piece(s) / pallet 420000
Depth of case 398 mm
Width of case 398 mm
Height of case 250 mm
Case volume 0.0396 cbm
Weight of product 0.0007 kg
Weight of case 7.36 kg
EAN of inner box 4038917084222
EAN case 4038917003087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

