



Product description

Order number 95.1007.001
Product description Tourniquet, material: PL, latex-free, 1 piece(s)/bag

Product characteristics

Colour of print/label white
Application venous blood collection

Material & colours

Product material polyester (PL)
Colour of product red

Purity & certification

Product category Medical device
CE certified CE supplier
Latex-free yes
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	200
Piece(s) / pallet	12800
Depth of case	330 mm
Width of case	330 mm
Height of case	220 mm
Case volume	0.024 cbm
Weight of product	0.0306 kg
Weight of case	6.86 kg
EAN of inner box	4260150884707
EAN case	4260150884806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

