



### **Product description**

Order number	14.1207
Product description	Blood Culture Adapter LongNeck, For narrow bottle neck, membrane screw cap assembled, connection: blood culture bottle, for use with (Safety) Multifly® needle, 1 piece(s)/blister, 50 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

### **Product characteristics**

Design	membrane screw cap assembled
Suitable for	use with (Safety) Multifly® needle
Connection	blood culture bottle

### **Size**

Diameter of opening	19 mm
Length of product	74 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Colour of product                      transparent

### **Purity & certification**

Product category                      In vitrodiagnostic device, CE  
CE certified                              CE  
Purity standard                         sterile, non-pyrogenic/endotoxin-free  
Sterilisation                             Electron irradiation  
Batched                                    yes

### **Packaging**

Minimum order qty.                    250  
Type of smallest                         blister  
subpackaging  
Piece(s) / inner box                    50  
Piece (s) / outer case                 250  
Piece(s) / pallet                         9000  
Depth of case                            578 mm  
Width of case                            368 mm  
Height of case                           195 mm  
Case volume                             0.0415 cbm  
Weight of product                      0.0041 kg  
Weight of case                           2.14 kg  
EAN of inner box                        4038917364188  
EAN case                                 4038917364171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

