



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 51.1638.940  |
| Product description | Safety-Multifly® BKF set 21G, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 200 mm, with blood culture adapter set, set components: 1 Safety-Multifly® needle 21G with long tubing, DEHP-free (REF 85.1638.235) and 1 Blood Culture Adapter Universal (REF 14.1207), 1 piece(s)/blister, sterile, 1 piece(s)/case |

### **Product characteristics**

|                |                           |
|----------------|---------------------------|
| Design         | blood culture adapter set |
| Coding colour  | green                     |
| Safety product | yes                       |
| Connection     | blood culture bottle      |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety-Multifly® BKF set 21G, 21G x 3/4", green,**  
**tube length: 200 mm, 1 piece(s)/blister**

**Page 2**

**Size**

|                               |               |
|-------------------------------|---------------|
| Diameter of needle (external) | 21 G / 0,8 mm |
| Diameter of tube              | 2.3 mm        |
| Length of needle              | 19 mm         |
| Length of tube                | 200 mm        |

**Purity & certification**

|                  |                |
|------------------|----------------|
| Product category | Medical device |
| Purity standard  | sterile        |
| Batched          | yes            |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 50000         |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 50            |
| Piece(s) / pallet             | 3200          |
| Depth of case                 | 398 mm        |
| Width of case                 | 298 mm        |
| Height of case                | 215 mm        |
| Case volume                   | 0.0255 cbm    |
| Weight of product             | 0.0082 kg     |
| Weight of case                | 1.09 kg       |
| EAN of inner box              | 4038917408738 |
| EAN case                      | 4038917408721 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**