



**Product Specification**  
**S-Monovette® needle, 21G x 1 1/2", green, 1  
piece(s)/blister**

**Page 1**



**Product description**

Order number	85.1162
Product description	S-Monovette® needle, needle (ØxL): 21G x 1 1/2", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

**Product characteristics**

Coding colour	green
Safety product	no
Connection	membrane screw cap

**Size**

Diameter of needle (external)	21 G / 0,8 mm
Length of needle	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® needle, 21G x 1 1/2", green, 1  
piece(s)/blister**

**Page 2**

**Purity & certification**

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0.0252 cbm
Weight of product	0.0017 kg
Weight of case	2.89 kg
EAN of inner box	4038917085526
EAN case	4038917004176

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**