



**Product Specification**  
**Multifly® needle, 21G x 3/4", green, tube length:  
80 mm, 1 piece(s)/blister**

**Page 1**

### Product description

Order number	85.1638.005
Product description	Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 80 mm, with Multi-adapter, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case

### Product characteristics

Design	Multi-adapter
Coding colour	green
Safety product	no
Connection	membrane screw cap

### Size

Diameter of needle (external)	21 G / 0,8 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	80 mm

### Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

### Packaging

Minimum order qty.	480
Type of smallest subpackaging	blister
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Multifly<sup>®</sup> needle, 21G x 3/4", green, tube length:  
80 mm, 1 piece(s)/blister**

Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.0027 kg
Weight of case	3.19 kg
EAN of inner box	4038917423762
EAN case	4038917423472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

