



Product Specification
Sample tube, Serum Gel CAT, 4.4 ml, cap brown,
(LxØ): 75 x 13 mm, with paper label



Product description

Order number	32.311.001
Product description	Sample tube, nominal volume: 4.4 ml, cap brown, preparation: Clotting activator/gel, push cap, (LxØ): 75 x 13 mm, round base, with paper label, 100 piece(s)/bag

Product characteristics

Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Closure type	push cap
Base shape	round base
Recommended centrifugation	2500 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, Serum Gel CAT, 4.4 ml, cap brown,
(LxØ): 75 x 13 mm, with paper label

Page 2

Size

Sample volume	4.4 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0043 kg
Weight of case	4.69 kg
EAN of inner box	4038917105262
EAN case	4038917010412

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT