



**Product Specification**  
**Sample tube, EDTA K3E, 5 ml, cap red, (LxØ): 75 x 13 mm, with print**

**Page 1**



**Product description**

Order number	32.332
Product description	Sample tube, nominal volume: 5 ml, cap red, preparation: K3 EDTA, push cap, (LxØ): 75 x 13 mm, round base, with print, 100 piece(s)/bag

**Product characteristics**

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	red
Graduation	yes
Cap	Skirted push cap
Closure type	push cap
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, EDTA K3E, 5 ml, cap red, (LxØ): 75 x 13 mm, with print**

**Page 2**

**Size**

Sample volume	5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0036 kg
Weight of case	4.12 kg
EAN of inner box	4038917081184
EAN case	4038917000321

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**