



### **Product description**

Order number	41.1392.005
Product description	Micro sample tube Serum CAT, venous blood collection, preparation: Clotting activator, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: white, colour code EU, with paper label, label/print: black, skirted conical base, 100 piece(s)/bag

### **Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with paper label
Colour of print/label	black
Graduation	yes
Closure type	screw cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,3 ml
Recommended centrifugation	10000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Centrifugation time 5 Min

### **Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 47 mm  
Length including cap 47 mm  
Length excluding cap 45 mm

### **Material & colours**

Colour of product transparent  
Closure material Polypropylene (PP)  
Colour of cap white  
Tube material Polypropylene (PP)

### **Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0016 kg  
Weight of case 1.93 kg  
EAN of inner box 4038917128261  
EAN case 4038917021623

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

