



Product Specification
Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print

Page 1



Product description

Order number 47.556
Product description Sample tube, nominal volume: 3 ml, cap red, preparation: K3 EDTA, screw cap, (LxØ): 82 x 11.5 mm, flat base, with print, 100 piece(s)/bag

Product characteristics

Type of preparation K3 EDTA
Preparation concentration 1.6 mg/ml blood
Colour code EU
Label/ Print with print
Colour of print/label red
Graduation yes
Closure type screw cap
Base shape flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Sample tube, EDTA K3E, 3 ml, cap red, (LxØ): 82 x 11.5 mm, with print

Size

Sample volume	3 ml
Diameter	11.5 mm
Length of product	82 mm
Length including cap	82 µm
Length excluding cap	66 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0027 kg
Weight of case	6.2 kg
EAN of inner box	4038917082440
EAN case	4038917001533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

