



Product Specification
Sample tube, EDTA K3E, 5 ml, cap red, (LxØ): 57 x 16.5 mm, with paper label

Page 1



Product description

Order number	49.355.001
Product description	Sample tube, nominal volume: 5 ml, cap red, preparation: K3 EDTA, screw cap, (LxØ): 57 x 16.5 mm, flat base, with paper label, 100 piece(s)/bag

Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Closure type	screw cap
Base shape	flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Sample tube, EDTA K3E, 5 ml, cap red, (LxØ): 57 x 16.5 mm, with paper label

Page 2

Size

Sample volume	5 ml
Diameter	16.5 mm
Length of product	57 mm
Length excluding cap	57 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0043 kg
Weight of case	4.8 kg
EAN of inner box	4038917105620
EAN case	4038917010726

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT