



Product description

Order number	01.1602.001
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 7.5 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with paper label

Page 2

Connection membrane screw cap

Size

Sample volume 7.5 ml
Outer diameter of cap 18 mm
Diameter 15 mm
Length including cap 108 mm
Length excluding cap 92 mm
Length including cap and plunger 129 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap brown
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 24000
Depth of box 189 mm
Width of box 134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,
(LxØ): 92 x 15 mm, with paper label**

Page 3

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0088 kg
Weight of case	5.24 kg
EAN of inner box	4038917131797
EAN case	4038917023979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT