



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 8.2 ml, cap green, (LxØ): 92 x 15 mm, with paper label

Page 1



Product description

Order number	01.1606.001
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 8.2 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 8.2 ml, cap green, (LxØ): 92 x 15 mm, with paper label

Liquid-tight Connection yes
 membrane screw cap

Size

Sample volume 8.2 ml
 Outer diameter of cap 18 mm
 Diameter 15 mm
 Length including cap 108 mm
 Length excluding cap 92 mm
 Length including cap and plunger 129 mm

Material & colours

Colour of product transparent
 Closure material High Density Polyethylene (HD-PE)
 Colour of cap green
 Membrane material rubber
 Tube material Polypropylene (PP)
 Piston material High Density Polyethylene (HD-PE)
 Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
 Product category In vitrodiagnostic device, CE
 CE certified CE
 Purity standard sterile
 Sterilisation Electron irradiation
 Batched yes

Packaging

Minimum order qty. 500
 Type of smallest subpackaging case
 Piece(s) / inner box 50
 Piece (s) / outer case 500
 Piece(s) / pallet 32000
 Depth of box 189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 8.2 ml, cap green, (LxØ): 92 x 15 mm, with paper label

Page 3

Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0085 kg
Weight of case	5.07 kg
EAN of inner box	4038917132046
EAN case	4038917024150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT